
Pharmamed
MEDICAL DEVICES

BILLROTH II SPHINCTEROTOME

S-Shape | Pull type



CODE	CUTTING WIRE	GUIDE WIRE	WC	PRODUCT
EST-02-T50-20	20 MM	0,035	2,8 MM	BILROTH II SPINCTERATOME 20



Billroth II sphincterotome

Used for cannulation of the ductal system and for sphincterotomy in patients with Billroth II anatomy.
It maintains its S-shape thanks to the memory wire inside. It automatically orients itself in the correct cutting direction with the duodenoscope, gastroscope or colonoscope.
Supplied sterile and is disposable. Intended for single use only.