

Pharmamed
MEDICAL DEVICES

ELECTROTOM SFINKTERATOM









Used for cannulation of the ductal system and for sphincterotomy. Sphincterotome is designed to aid navigation through the papilla and the fronds, potentially reducing some risks associated with sphincterotomy. And, with the integrated hub, you have the ability to inject while the wire is in the lumen, giving you the efficiency and procedural control you've come to expect from the device. Supplied sterile and is disposable. Intended for single use only.

Triple lumen sphincterotome. Separate lumens for guidewire, cutting wire and injection of contrast medium.
Dome tip. Easy access, orientation and security.
Ergonomic, rotating handle. Engineered to aid precise tip movement and wire manipulation to facilitate cannulation.
Improved guidewire port. You can easily advance even a J-shaped guidewire without feeling any snagging.
Markers printed with PTFE ink. Further reduces friction for easier cannulation.
Unique stopper for needle knife sphincterotome. Strictly limits the blade output (Default length 5mm).

CODE	CUTTING WIRE	GUIDE WIRE	WC	PRODUCT
EST-01-T50-20	20 MM	0,035	3,2 MM	SPINCTERATOME 20
EST-01-T50-25	25 MM	0,035	3,2 MM	SPINCTERATOME 25
EST-01-T50-30	30 MM	0,035	3,2 MM	SPINCTERATOME 30
EST-01-T50- I30	30 MM	0,035	3,2 MM	SPINCTERATOME 30 Partially insulated cutting wire